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## Determination of Substances of Concern Regulated in EU

Sumika Chemical Analysis Service, Ltd.  
Ehime Laboratory

Hiromu OGASAWARA  
Kei TANAKA  
Shuichiro MANABE

In recent years, advanced regulations concerning restrictions on the use of hazardous substances have been enforced in the EU, as environmental pollution control acts. Subsequent to the ELV for automobiles, restrictions on the use of hazardous substances are being examined for the field of electric and electronic equipment, where both WEEE and RoHS directives have been enacted. In the ELV, the use of four heavy metals (Cadmium, Lead, Mercury, Chromium) has been restricted. The use of the same heavy metals as in the ELV, as well as two additional flame retardants (Polybrominated biphenyl, Polybrominated diphenyl ether) has been restricted in RoHS. Under these circumstances, the control of certain hazardous substances in products takes on great importance.

This paper describes the use of various analytical equipment in screening analysis by XRF and precision analysis for Cd, Pb, Hg, Cr, Br.

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Table 1 Restrictions on substances of concern

country/area	
EU	Waste electrical and electronic equipment (WEEE) Restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) End of life vehicles (ELV)
Japan	The marking for presence of the specific chemical substances for electrical and electronic equipment (J-Moss)
USA (California)	Electronic Waste Recycling Act (S.B.20) Electronic Waste, Advanced Disposal Fees (S.B.50)
USA (CONEG)	The Model Toxics in Packaging Legislation <sup>1)</sup>
Korea	The Act for Resource Recycling of Electrical/Electronic Products and Automobiles (Korean WEEE, RoHS)
China	Management Method on the Prevention and Control of Pollution Caused by Electronic Information Products (Chinese RoHS)

**Table 2** Maximum concentration values of substances of concern


Substance	ELV	RoHS	unit : mg/kg	
			Procurement guideline	Procurement guideline
Heavy metal	Pb	1000	100	250
	Hg	1000	—	1000
	Cd	100	5	50
	Cr <sup>6+</sup>	1000	—	1000
Flame retardant	PBB	1000	—	1000
	PBDE	1000	—	1000

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**Table 3** Screening limits for regulated elements

Element	Lower limit	Upper limit
Cd, Pb, Hg	$C \times (1 - S) - n\sigma$	$C \times (1 + S) + n\sigma$
Br, Cr	$C \times (1 - S) - n\sigma$	—

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
Fig. 3 Calibration of elements in standard sample of PVC

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
Fig. 4 Comparison of X-ray strength in different matrix

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**Table 5** Analytical results of Cd,Cr,Hg,Pb and Br in polyethylene (BCR-680)

Element	Certified value	unit : mg/kg				
		ICP-MS	ICP-AES	CV-AAS	ET-AAS	IC
Cd	141	145	143	—	—	—
Cr	115	122	110	—	—	—
Hg	25	26	—	25	—	—
Pb	108	109	—	—	113	—
Br	808	—	—	—	—	792

Analysis : Pretreatment

ICP-MS : Mountable acid digestion

ICP-AES : Demountable acid digestion

CV-AAS : Acid digestion attached reflux condenser

ET-AAS : Mountable acid digestion

IC : Combustion in the atmosphere of oxygen

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**Table 6** The classification of flame retardant


Type of compound	Compound
Inorganic	Antimony (III) Oxide, Magnesium hydroxide, Zinc borate
Bromide	Deca brominated biphenyl ether, Tetrabromobis phenol-A, Hexabromo benzene
Chloride	Chlorinated paraffine, Chlorinated polyethylene
Phosphate	Triphenyl phosphate, Tricresyl phosphate
Others	Triazine compounds, Guanidine compounds, Silicone oil

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**Fig. 6** The chemical formula of PBDE and PBB

**Table 7** Name of homologue and number of isomer of PBBs and PBDEs

Number of bromine	Name of homologue	Number of isomer
1	Mono brominated	3
2	Di brominated	12
3	Tri brominated	24
4	Tetra brominated	36
5	Penta brominated	42
6	Hexa brominated	36
7	Hepta brominated	24
8	Octa brominated	12
9	Nona brominated	3
10	Deca brominated	1
	Total	209

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